

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\\Research\\Rorden_Prisma\\MEMPRAGE\\Localizer
 TA:0:22 PAT:Off Voxel size:0.5×0.5×7.0 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	On
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	3
Slices	3
Dist. factor	300 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	7.0 mm
TR	8.6 ms
TE	4.00 ms
Averages	2
Concatenations	5
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	91 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
TD	0 ms
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Nr. of slice groups	3
Slices	3
Dist. factor	300 %
Position	L0.0 A30.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.255798 MHz
Correction factor	1
SRFExcit 1H	42.652 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	2D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	320 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Liver registration	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPIMAGE\T1_mprage_ns_sag_p2
 TA:6:17 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A14.0 F34.8 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2250.0 ms
TE	4.11 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	925 ms
Flip angle	9 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	80
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A14.0 F34.8 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
Multi-slice mode	Single shot
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE1	Off
HE3	Off
NE1	Off
HE2	Off
HE4	Off
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 A14.0 F34.8 mm
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	192 mm

Frequency 1H	123.255798 MHz
Correction factor	1
SLoopIRns1 1H	403.002 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	925 ms
Dark blood	Off
Resp. control	Off
Inline	
Inline Composing	Off
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Bandwidth	150 Hz/Px
Flow comp.	No
Echo spacing	9.3 ms
Turbo factor	192
RF pulse type	Normal
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPORAGE\tfl_mgh_multiecho
 TA:6:03 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl_me

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2530.0 ms
TE 1	1.61 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1100 ms
Flip angle	7.0 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	176
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE2	Off
HE4	Off
NE2	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	176 mm
Frequency 1H	123.255798 MHz
Correction factor	1
SLoopIRns1 1H	403.002 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	4.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1100 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Bandwidth 1	650 Hz/Px
Flow comp. 1	No
Echo spacing	9.5 ms
Turbo factor	176
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Gradient spoiling	Siemens
Gradient moment factor	1
Averaging	None
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	4
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPIMAGE\tfl_mgh_multiecho_fat
 TA:6:03 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl_me

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2530.0 ms
TE 1	1.61 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1150 ms
Flip angle	7.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	176
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE2	Off
HE4	Off
NE2	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	176 mm
Frequency 1H	123.255798 MHz
Correction factor	1
ExcitWEs 0 1H	52.306 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	4.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1150 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Bandwidth 1	650 Hz/Px
Flow comp. 1	No
Echo spacing	10.7 ms
Turbo factor	176
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Gradient spoiling	Siemens
Gradient moment factor	1
Averaging	None
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	4
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPIMAGE\tfl_mgh_multiecho_fat_HCP
 TA:7:24 PAT:2 Voxel size:0.8×0.8×0.8 mm Rel. SNR:1.00 :tfl_me

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	0.80 mm
TR	2530.0 ms
TE 1	1.36 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1150 ms
Flip angle	8.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE2	Off
HE4	Off
NE2	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	167 mm
Frequency 1H	123.255798 MHz
Correction factor	1
ExcitWEs 0 1H	59.778 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	4.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1150 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Contrasts	4
Bandwidth 1	650 Hz/Px
Flow comp. 1	No
Echo spacing	10.8 ms
Turbo factor	208
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Gradient spoiling	Siemens
Gradient moment factor	1
Averaging	None
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	4
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPAGENT1_HCP
 TA:6:38 PAT:2 Voxel size:0.8×0.8×0.8 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A11.5 F0.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	23.1 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	2400.0 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1060 ms
Flip angle	8 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A11.5 F0.6 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	23.1 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 A11.5 F0.6 mm
Rotation	0.00 deg
F >> H	256 mm
A >> P	240 mm
R >> L	167 mm
Frequency 1H	123.255798 MHz
Correction factor	1
ExcitWEns 0 1H	59.778 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1060 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	210 Hz/Px
Flow comp.	No
Echo spacing	8.1 ms
Turbo factor	256
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\\Research\\Rorden_Prisma\\MEMPRAGE\\T1_mprage_ns_sag_p2_fat
 TA:6:17 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A14.0 F34.8 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2250.0 ms
TE	4.11 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1040 ms
Flip angle	9 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	80
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A14.0 F34.8 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
Multi-slice mode	Single shot
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE1	Off
HE3	Off
NE1	Off
HE2	Off
HE4	Off
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 A14.0 F34.8 mm
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	192 mm

Frequency 1H	123.255798 MHz
Correction factor	1
ExcitWEns 0 1H	67.250 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1040 ms
Dark blood	Off
Resp. control	Off
Inline	
Inline Composing	Off
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Bandwidth	150 Hz/Px
Flow comp.	No
Echo spacing	10.6 ms
Turbo factor	192
RF pulse type	Normal
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPRAGE\Localizer
 TA:0:22 PAT:Off Voxel size:0.5×0.5×7.0 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	On
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	3
Slices	3
Dist. factor	300 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	7.0 mm
TR	8.6 ms
TE	4.00 ms
Averages	2
Concatenations	5
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	91 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
TD	0 ms
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Nr. of slice groups	3
Slices	3
Dist. factor	300 %
Position	L0.0 A30.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.255798 MHz
Correction factor	1
SRFExcit 1H	42.652 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	2D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	320 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Liver registration	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPRAGE\me_FieldMap_GRE
TA:1:46 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :fm_r

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	36
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	400.0 ms
TE 1	5.19 ms
Averages	1
Concatenations	2
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magn./Phase
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	36
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	129 mm
Frequency 1H	123.255798 MHz
Correction factor	1
01GreFCE 1H	127.955 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	2D
Averaging mode	Long term
Multi-slice mode	Interleaved
Asymmetric echo	Off
Contrasts	2
Bandwidth	260 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Distortion Corr.	Off
Contrasts	2

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPRAGE\me_fmRI
 TA:0:20 PAT:2 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	30
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.00 mm
TR	2000 ms
TE 1	20.00 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	79 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	6
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Single-shot
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	30
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	108 mm
Frequency 1H	123.255798 MHz
Correction factor	1
SincExcRF 1H	92.995 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Contrasts	2
Bandwidth	2298 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.52 ms
EPI factor	64
Gradient mode	Performance

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	20
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off
Contrasts	2

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPRAGE\me_MPRAGE
 TA:6:03 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl_me

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2530.0 ms
TE 1	1.61 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1100 ms
Flip angle	7.0 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	176
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE2	Off
HE4	Off
NE2	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	176 mm
Frequency 1H	123.255798 MHz
Correction factor	1
SLoopIRns1 1H	403.002 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	4.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1100 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Bandwidth 1	650 Hz/Px
Flow comp. 1	No
Echo spacing	9.5 ms
Turbo factor	176
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Gradient spoiling	Siemens
Gradient moment factor	1
Averaging	None
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	4
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\MEMPIMAGE\me_MPRAGE_FAT
 TA:6:03 PAT:2 Voxel size:1.0×1.0×1.0 mm Rel. SNR:1.00 :tfl_me

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2530.0 ms
TE 1	1.61 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1150 ms
Flip angle	7.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	176
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
HE2	Off
HE4	Off
NE2	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	176 mm
Frequency 1H	123.255798 MHz
Correction factor	1
ExcitWEs 0 1H	52.306 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	4.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1150 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Contrasts	4
Bandwidth 1	650 Hz/Px
Flow comp. 1	No
Echo spacing	10.7 ms
Turbo factor	176
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Gradient spoiling	Siemens
Gradient moment factor	1
Averaging	None
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	4
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D**Table of contents**

```

\USER
| Research
|   | Rorden_Prisma
|   |   | MEMPRAGE
|   |   |   | Localizer
|   |   |   | T1_mprage_ns_sag_p2
|   |   |   | tfl_mgh_multiecho
|   |   |   | tfl_mgh_multiecho_fat
|   |   |   | tfl_mgh_multiecho_fat_HCP
|   |   |   | T1_HCP
|   |   |   | T1_mprage_ns_sag_p2_fat
|   |   |   | Localizer
|   |   |   | me_FieldMap_GRE
|   |   |   | me_fmRI
|   |   |   | me_MPRAGE
|   |   |   | me_MPRAGE_FAT

```